ALBERTA PRECISION LABORATORIES

Leaders in Laboratory Medicine

INITIAL SET OF LABORATORY TESTS in a patient with suspected viral haemorrhagic fever - EDMONTON ZONE ONLY

NO SAMPLE COLLECTION SHOULD OCCUR UNTIL APPROVED BY THE MEDICAL OFFICER OF HEALTH (MOH)

	Tests available ^a	Container and Required Volume(s) ^{b,c}
EVD testing	Ebola virus serology (NML)	1 serum separator tube (SST) (GOLD top) tube (3 ml)
	Ebola virus RT PCR* (ProvLab/NML)	2 EDTA (PURPLE top) tubes each with 3-4 ml
Malaria	Malaria rapid detection (ProvLab UAH) Malaria PCR (ProvLab UAH)	1 EDTA (PURPLE top) tube with 3-5 ml
Hematology	CBC with differential (ProvLab UAH): Haemoglobin, Hematocrit, platelets, red blood cell count and white blood cell count with 3 part differential (neutrophils, lymphocytes and other)	2 EDTA (PURPLE top) 4 ml tube • Required <u>each</u> time testing is needed.
Chemistry	General Chemistry (ProvLab UAH): Creatinine and urea, electrolytes (Na+, K+, Cl, HCO ₃ -), ionized Calcium, glucose, blood gases, lactate, calcium, liver enzymes (ALT, ALP, AST, GGT), total bilirubin, amylase, total protein Serum Beta-HCG, Troponin I	 2 Lithium Heparin (GREEN top) fill to capacity Required <u>each</u> time testing is needed.
Blood cultures	Adult and pediatric blood cultures (ProvLab UAH)	 Adult Blood cultures (Bactec ☑): 2 sets (Minimum 3 bottles: 2 aerobic and one anaerobic bottle), each with 10 ml Pediatric Blood cultures (Bactec): ≥ 14 years old: 1 pediatric and 1 anaerobic bottle < 14 years old: 1 pediatric bottle

^a All testing done at ProvLab UAH or sent from ProvLab to the National Microbiology Laboratory.

^b Ebola RT PCR testing performed within 3 days of fever onset may not have detectable Ebola virus. Repeat Ebola PCR testing on a new blood specimen ≥3 days after onset of fever to confirm negative results. ^cThis is the <u>ONLY</u> testing available for a suspect/confirmed Ebola patient. No other testing is available.

Questions/concerns for lab testing/results for VHF patients:

• Virologist-on-call at the Provincial Laboratory: page via 780-407-8200